

FLINN SCIENTIFIC INC.

"Your Safer Source for Science Supplies"

Material Safety Data Sheet (MSDS)

MSDS #: 528.00

Revision Date: November 25, 2002

Section 1 — Chemical Product and Company Identification

Millon's Reagent Solution

Flinn Scientific, Inc. P.O. Box 219 Batavia, IL 60510 (800) 452-1261

CHEMTREC Emergency Phone Number: (800) 424-9300

Section 2 — Composition, Information on Ingredients

Mercuric Nitrate (10415-75-5) 38%, Nitric Acid (7697-370-2) 18%, and Water (7732-18-5) 44%.

CAS#: None Established

Section 3 — Hazards Identification

Colorless liquid with characteristic nitric acid odor.

Corrosive liquid, severe body tissue hazard. Avoid contact with any body tissues. Highly toxic as mercury. Irritating to respiratory system. Avoid inhalation.

FLINN AT-A-GLANCE

Health-3
Flammability-0
Reactivity-2
Exposure-3
Storage-3

0 is low hazard, 3 is high hazard

Section 4 — First Aid Measures

Call a physician, seek medical attention for further treatment, observation and support after first aid.

Inhalation: Remove to fresh air at once. If breathing has stopped give artificial respiration immediately.

Eye: Immediately flush with fresh water for 15 minutes.

External: Wash continuously with fresh water for 15 minutes.

Internal: Give large quantities of water. Call a physician or poison control at once.

Section 5 — Fire Fighting Measures

Not combustible liquid.

When heated to decomposition, emits toxic fumes of mercury and NOx.

Fire Fighting Instructions: Use triclass, dry chemical fire extinguisher. Firefighters should wear PPE and SCBA with full facepiece operated in positive pressure mode.

NFPA CODE

None Established

Section 6 — Accidental Release Measures

Restrict unprotected personnel from area. Remove all ignition sources and ventilate area. Contain spill with sand and absorbent material, neutralize with sodium bicarbonate or calcium hydroxide and deposit in sealed bag or container. See Sections 8 and 13 for further information.

Section 7 — Handling and Storage

Flinn Suggested Chemical Storage Pattern: Inorganic #9. Store with inorganic acids.

Store in a dedicated acid cabinet and away from any source of water; if an acid cabinet is not available, store in Flinn Saf-Cube.

Store in a cool dry place. Use and dispense in a hood.

Section 8 — Exposure Controls , Personal Protection

Avoid contact with eyes, skin and clothing. Wear chemical splash goggles, chemical-resistant gloves and chemical-resistant apron.

Use exhaust ventilation to keep airborne concentrations low.

Section 9 — Physical and Chemical Properties

Colorless liquid with characteristic nitric acid odor.

Section 10 — Stability and Reactivity

Avoid contact with ethyl alcohol, sulfur and exposure to light.

Shelf life: Good, if stored safely.

Section 11 — Toxicological Information

Acute effects: Highly toxic, corrosive

Chronic effects: N.A.

Target organs: N.A.

ORL-RAT LD50: 297 mg/kg as mercury

IHL-RAT LC50: N.A.

SKN-RBT LD50: N.A.

N.A. = Not available, not all health aspects of this substance have been fully investigated.

Section 12 — Ecological Information

Data not yet available.

Section 13 — Disposal Considerations

Please consult with state and local regulations.

Flinn Suggested Disposal Method #27f is one option.

Section 14 — Transport Information

Shipping Name: Not regulated

Hazard Class: N/A

UN Number: N/A

N/A = Not applicable

Section 15 — Regulatory Information

Not listed.

Section 16 — Other Information

Consult your copy of the Flinn Scientific Catalog/Reference Manual for additional information about laboratory chemicals. This Material Safety Data Sheet (MSDS) is for guidance and is based upon information and tests believed to be reliable. Flinn Scientific Inc. makes no guarantee of the accuracy or completeness of the data and shall not be liable for any damages relating thereto. The data is offered solely for your consideration, investigation, and verification. Flinn Scientific Inc. assumes no legal responsibility for use or reliance upon this data.

FLINN SCIENTIFIC INC.

"Your Safer Source for Science Supplies"

Flinn Chemicals--Your Best Choice

flinn@flinnsci.com www.flinnsci.com
P.O. Box 219 Batavia IL 60510
(800) 452-1261 Fax (866) 452-1436